

# The TF20 microscope interfaces

## The microscope console

**Location** : right in front of you on the desk (screen on the left)

**Role:** TEM User Interface provides direct control of some of the most important microscope parts:

- Magnification
- Electron dose: spot size (C1 lens) and beam spread (C2 lens)
- Focus (objective lens current)
- Stage position (with the joystick and Z buttons)
- Beam position (with the trackball)
- and many other actions through the programmable switches (L/R1-2-3) and Multifunction X/Y knobs.



\* [For more details see:](#)

\*

### TUI interface

TUI interface is divided in 4 zones

- Zone1 : TUI tabbed pages hosting different commands allowing to control the ME and to know its working state.
- Zone 2 : Window giving values of working state : HT value, stage position, magnification, etc ...
- Zone 3 : Additional TUI Tabbed page
- Zone4 : Window defining functions of control panel button:



## Gatan DigitalMicrograph User Interface (DM)

**Location:** On the left screen

**Role:** It controls the K2 direct electron detector.

- The beginner user will not use it.
- [The expert user will use it to:](#)
  - prepare Dark, linear and counting / Super Resolution gain references.

- measure the dose.
  - monitor the K2 health status.
  - troubleshoot the system.
  - perform cryo-cycle of the detector
- 

## SerialEM User Interface (SEM)

**Location:** On the left screen

**Role:** It provides camera/microscope commands to allow the user to screen its grids and setup an automated data collection.

- **The beginner user will use it to:**
    - record the Atlases.
    - select squares and prepare the square maps at eucentric height.
    - select targets (holes).
    - screen targets (holes).
  - **The expert user will use it to:**
    - perform all the actions described above.
    - perform beam corrections (astigmatism, coma-free).
    - perform calibrations.
    - setup the dose rate, total dose and dose fractionning.
    - setup the time and imaging modes for the different beam modes.
    - center the different beam modes all together on the camera axis.
- 

## The TEM User Interface (TEM)

**Location:** On the middle screen

**Role:** It controls all the microscope systems.

- **The beginner user will use it to:**
    - dock/undock the cassette from the NanoCab to the autoloader.
    - load/unload the grids in the column.
    - change and center the Condenser 2 and objective apertures.
    - control the autoloader turbo pump.
    - perform beam-shift adjustments.
    - manually set eucentric height.
  - **The expert user will use it to:**
    - perform all the actions described above.
    - align the microscope.
-

- troubleshoot the system.
- manage the cold system.
- perform cryo-cycle of the column/autoloader.

---

## The FluCam User Interface (FC)

**Location:** On the right screen

**Role:** It allows the visualization of the electron beam and the sample image.

- **The beginner user will use it to:**
  - center the beams.
  - center Condenser 2 and objective apertures.
  - look for and move to particular features on their grids.
  - manually set eucentric height.
- **The expert user will use it to:**
  - perform all the actions described above.
  - align the microscope.

From:

<https://bsi.inscog.eu/> - **BSI wiki**

Permanent link:

[https://bsi.inscog.eu/doku.php?id=glacios\\_interface&rev=1679328309](https://bsi.inscog.eu/doku.php?id=glacios_interface&rev=1679328309)

Last update: **2023/11/01 20:15**

