

Expression & purification of 6his-SmPRMT3(1-553)

Construct

MGSSHHHHHHSSGTGSGENLYFQGHMSPCDQSSNYDLDMEPNDSSYFGSYGHFEIHGEMINDRVRTESYVNFILSN
AEKYFKHKIILDVGSVSGSILSIIAAQAGASHVYGVEAADEIYAASHETLRVNNLLERVTFIHGQAESVELPVKKVDVIISE
WMGYFLFFESMLDSVLKMASKYLSRDGHIFPRHYTLNLLGVQCSEQLRKRRLEHWNNVYGYNMPALRRAALSEAHV
LNLTNEHVTPPISPITILTQSFELVALDLDDMHRNRIYNLSNHCSLLCEQKFHLTIQPTTDINNNSSSSSYELDAIVGYFD
VRFDDADCKVEFSTSPPTPLTHWKQTLFLDKPIRVKPGDKISGIITIRRATTDNRGLEINLLIGETENSLEIKQTFDLIG

Length = 396 AA

Molecular weight = 44782.3 Da

pI 5.66

$\xi(\text{red})$ 40340 L.M⁻¹

Expression

Grow pnEAvH_SmPRMT3-184 BL21(DE3) transformants on an LB-Amp plate (1 plate/L culture), overnight 37°C. Resuspend cells in LB medium (5 mL/plate) and inoculate liquid **2X** LB medium + 100 µg/mL⁻¹ Ampicilline. Grow culture(s) at 37°C, 200 rpm and measure OD600 every 20 minutes. When OD600 reaches 0.4, set temperature to 20°C. When OD600 reaches 0.6 to 0.8, induce 6his-SmPRMT3(184-553) with 1 mM IPTG and leave the cultures to express the protein overnight.

Purification (for 3L)

Pellet cells by centrifugation (4000 rpm, 45 minutes), then resuspend pellets in 200 mL Lysis buffer

Buffers

Lysis buffer: 50 mM BTP pH9, 250 mM NaCl, 1 mM TCEP, 0.01% NP40, Complete[®]-EDTA free protease inhibitor cocktail (1 pill / 50 mL)

From:
<https://bsi.inscog.eu/> - BSI wiki

Permanent link:
https://bsi.inscog.eu/doku.php?id=purification:expression_and_purification_of_smprmt3_184-553&rev=1477662145

Last update: 2023/11/01 20:17

