

Expression & purification of 6his-SmPRMT3(1-553)

Construct

MGSSHHHHHHSSGTGSGENLYFQGHMSPCDQSSNYDLDMEPNDSSYFGSYGHFEIHGEMINDRVRTESYVNFILSN
AEKYFKHKIILDVGS GSGILSIIAAQAGASHVYGVEAADEIYAASHETLRVNNLLERVTFIHGQAESVELPVKKVDVIISE
WMGYFLFFESMLDSVLKMASKYLSRDGHIFPRHYTLNLLGVQCSEQLRKRRLHWNVYGYNMPALRRAALSEAHV
LNLNEHVTPPISPITILTQSFELVALDLDMMHRNRINYLSNHCSLLCEQKFHLTIQPTTDINNNSSSSSYELDAIVGYFD
VRFDDADCKVEFSTSPPTPLTHWKQTLFLDKPIRVKPGDKISGIITIRRATTDNRGLEINLLIGETENSLEIKQTFDLIG

Length = 396 AA

Molecular weight = 44782.3 Da

pI 5.66

$\xi(\text{red})$ 40340 L.M⁻¹

Expresion

Grow pnEAvH_SmPRMT3-184 BL21(DE3) transformants on an LB-Amp plate (1 plate/L culture), overnight 37°C. Resuspend cells in LB medium (5 mL/plate) and inoculate liquide **2X** LB medium + 100 µg/mL⁻¹ Ampicilline. Grow culture(s) at 37°C, 200 rpm and measure OD600 every 20 minutes. When OD600 reaches 0.4, set temperature to 20°C. When OD600 reaches 0.6 to 0.8, induce 6his-SmPRMT3(184-553) with 1 mM IPTG and leave the cultures to express the protein overnight.

Purification (for 3L)

Pellet cells by centrifugation (4000 rpm, 45 minutes), then resuspend pellets in 200 mL Lysis buffer

Buffers

Lysis buffer (200 mL): 50 mM BTP pH9, 250 mM NaCl, 1 mM TCEP, 0.01% NP40, Complete® -EDTA free protease inhibitor cocktail (1 pill / 50 mL)

Talon wash buffer (100 mL): 50 mM BTP pH9, 250 mM NaCl, 1 mM TCEP, 10 mM imidazole

Talon elution buffer (20 mL): 50 mM BTP pH9, 250 mM NaCl, 1 mM TCEP, 150 mM imidazole

Dilution buffer (180 mL): 20 mM BTP pH9, 1 mM TCEP, 1 mM EDTA

ANX buffer (200 mL): 20 mM BTP pH9, 25 mM NaCl, 1 mM TCEP, 1 mM EDTA

Gradient buffer (100 mL): 20 mM BTP pH9, 1 M NaCl, 1 mM TCEP, 1 mM EDTA

GF buffer (350 mL): 20 mM BTP pH9, 50 mM NaCl, 1 mM TCEP, 1 mM EDTA

From:
<https://bsi.inscog.eu/> - BSI wiki

Permanent link:
https://bsi.inscog.eu/doku.php?id=purification:expression_and_purification_of_smprmt3_184-553&rev=1477662529

Last update: 2023/11/01 20:17

