

If you have any data coming from the microscopes, it should be available to transfer from
/mnt/gateway/teams/<your_team>/<your_project>/<your_dataset>

The best way to transfer it is to use rclone, here is a usage exemple:

```
rclone copy -P -l --transfers=6 /mnt/gateway/teams/<source>/  
/mnt/storage/teams/<destination>/
```

rclone behaves exactly like rsync but adds the parallel transfer feature `--transfers=6` being the number of parallel transfers, you can try to tweak it to optimize the transfer speed but 6 seems to work well.

From:

<https://bsi.inscog.eu/> - **BSI wiki**

Permanent link:

https://bsi.inscog.eu/doku.php?id=collected_datasets&rev=1683186625

Last update: **2023/11/01 20:14**

