

In-house TEV production

Clones

[pRK792 plasmid](#) - His-TEV(S219V)-Arg - Waugh

Sequence

GHHHHHHHGESLFGPRDYNPISSTICHLTNESDGHTTSLYGIGFGPFIITNKHLFRRNNGTLLVQSLHGVFKVKNTT
TLQQLIDGRDMIIRMPKDFPPFPQKLKFRPQREERICLVTTNFQTKSMSSMVSDTSCTFPSSDGIFWKHWIQTKDG
QCGSPLVSTRDGFIVGIHSASNFTNTNNYFTSVPKNFMELLTNQEAQQWVSGWRLNADSVLWGGHKVFMVKPEEPF
QPVKEATQLMNR**RRRRR**

His-TEV(L35V/S135G/S219V)-Arg - Cabrita¹⁾

His-TEV plasmid - Birck

pBS2440 optimized TEV plasmid - Seraphin

Protocols

[Waugh protocol](#)

[Cabrita protocol](#)

[Birck protocol](#)

[Seraphin protocol](#)

[TEV FAQ](#)

¹⁾

[Cabrita et al. \(2007\) Protein Science 16, 2360-2367](#)

From:

<https://bsi.inscog.eu/> - **BSI wiki**

Permanent link:

<https://bsi.inscog.eu/doku.php?id=purification:tev&rev=1468830322>

Last update: **2023/11/01 20:17**

